

Product Details

This product is a mixture of ten highly purified pre-stained proteins ranging from 8 kDa to 180 kDa. The 72 kDa band is orange, 8 kDa band is green, and other bands are blue. It is designed for observing protein separation during SDS-PAGE, verifying western transfer efficiency on membranes, and approximating the size of proteins.

Product Features

- Three Pre-Stained Color
- Neat, Bright and Well-distributed Bands
- More Accurate Molecular Weight Position

Storage Buffer

20 mM Tris-H3PO4 (pH 7.5), 2 mM EDTA, 1.5 % (W/V) SDS, 3 mM DTT, 0.1% (V/V) Proclin300, 15 % (V/V) Glycerol.

Shipping and Storage

The product is shipped with blue ice. Upon receipt, store it immediately at -20°C for long term storage.

This product is stable after storage at:

- -20°C for up to three years or 4°C for up to two months.

Procedure

1. Thaw the product at room temperature for a few minutes to dissolve precipitated solids. Do not boil!
2. Mix gently, but thoroughly, to ensure the solution is homogeneous.
3. Load the following volumes of the product on an SDS-PAGE:
 - 3-5 µL per well for mini gel
 - 5-10 µL per well for large gel

Use the same volumes for Western blot. The loading volumes listed above are recommended for gels with a thickness of 0.75-1.0 mm.

The loading volume should be doubled for 1.5 mm thick gels.


Product Specification

Number of Bands	10
Size Range	8 to 180 kDa
Stain Type	3 colors: Blue, Orange, Green
Molecular Weight	180, 130, 100, 72, 55, 43, 33, 25, 17, 8 kDa
Quantity	2 x 250 ul, 10 x 250 ul
System Type	SDS-PAGE, Western Blot

Notes

1. This product has been prepared in 1× SDS-PAGE loading buffer and can be used directly without boiling, diluting and adding reducing agent.
2. Longer transfer times or higher transfer voltages may be required for Western blot of large (>100 kDa) proteins.
3. Don't add SDS to transfer buffer. If SDS must be used, the concentration should not exceed 0.02-0.04%.
4. In low-percentage gels (< 10 %), the low-molecular weight proteins in the ladder may migrate with the dye front.
5. Pre-stained proteins can have different mobilities in various SDS-PAGE-buffer systems. However, they are suitable for approximate molecular weight determination when calibrated against unstained standards in the same system. See the table provided for migration patterns in different electrophoresis conditions.
6. For your safety and health, please wear a lab coat and disposable gloves.

SDS-PAGE



SDS-PAGE band profile of the Star Ribbon Pre-stained Protein Marker

Migration patterns of the Marker in different electrophoretic conditions

Gel type	Tris-Glycine						Bis-Tris						Tris-Acetate		Hepes-Tris	
	8%	10%	12.5%	15%	B4-20%	W4-20%	G4-12%	G8-16%	G4-20%	G4-12%	G8-16%	G4-20%	G10%	6%	T3-8%	W4-20%
Running buffer	Tris-Glycine						MES			MOPS			Tris-Acetate		Hepes	
Apparent Molecular Weights, kDa																
10	180	180	180	180	180	180	180	180	180	180	180	180	180	180	180	180
20	180	130	100	75	55	180	130	130	130	130	130	130	130	130	130	130
30	130	72	55	43	130	180	65	65	65	65	65	65	65	65	65	65
40	100	55	43	33	100	130	43	43	43	43	43	43	43	43	43	43
50	70	43	33	25	70	33	33	33	33	33	33	33	33	33	33	33
60	55	33	25	17	55	25	25	25	25	25	25	25	25	25	25	25
70	43	25	17	17	43	17	17	17	17	17	17	17	17	17	17	17
80	33	17	17	17	33	17	17	17	17	17	17	17	17	17	17	17
90	33	17	8	8	33	17	17	17	17	17	17	17	17	17	17	17
100	33	17	8	8	33	17	17	17	17	17	17	17	17	17	17	17

The apparent molecular weight of each protein (kDa) has been determined by calibration of each protein against an unstained protein ladder in specific electrophoresis conditions. Migration patterns were determined using commercial precast mini gels.