



## ELISA Sample Dilution Buffer -4

### Product information

Catalog Number: CEB-C302

Size: 15 mL/bottle

ELISA Sample Dilution Buffer-4 is a solution used to dilute and prepare samples for testing, facilitating the detection of the target substance with ClinMax™ ELISA reagents. It does not directly participate in the detection process itself.

### Storage

Keep the reagents stored at 2-8 °C. For opened kit and reconstituted reagents, can be stored for up to 30 days at 2-8 °C.

### NOTE:

This dilution buffer product is for research use only and is not for usage in diagnostic or therapeutic applications.

Please do not use the reagents after the expiration date indicated on the label.

Manufactured and distributed by  
ACRODiagnostics Inc.  
TEL: 010-67855298-8147

Distributed by:  
ACROBiosystems Inc.  
US & Canada TEL: +1 800-810-0816  
Asia & Pacific TEL: +86 400-682-252