



Synonym

capsid protein

Source

AAV2 VP1, Recombinant Protein(VP1-A5143) is expressed from E. coli cells. It contains AA Met1 - Leu735 (Accession # [P03135-1](#)).

Predicted N-terminus: Met

Molecular Characterization

Poly-his VP1(Met1 - Leu735)
P03135-1

This protein carries a polyhistidine tag at the N-terminus.

The protein has a calculated MW of 84.0 kDa. The protein migrates as 45-90 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE).

Endotoxin

Less than 1.0 EU per µg by the LAL method.

Purity

>85% as determined by SDS-PAGE.

Formulation

Supplied as 0.2 µm filtered solution in 20 mM Tris, 0.5 M NaCl, 0.5 M Arginine, pH8.0 with trehalose as protectant.

Contact us for customized product form or formulation.

Shipping

This product is supplied and shipped with dry ice, please inquire the shipping cost.


Storage

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- The product MUST be stored at -70°C or lower upon receipt;
- -70°C for 3 months under sterile conditions.

SDS-PAGE




AAV2 VP1, Recombinant Protein on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 85% (With [Star Ribbon Pre-stained Protein Marker](#)).

Bioactivity-ELISA

Discounts, Gifts,
and more!





Immobilized AAV2 VP1, Recombinant Protein (Cat. No. VP1-A5143) at 1 µg/mL (100 µL/well) can bind Anti-AAV,VP1+VP2+VP3,B1 Mouse Monoclonal,Purified,IgG1 with a linear range of 1-31 ng/mL (QC tested).

Background

Adeno-associated virus is a single-stranded DNA virus and the current scientific consensus is that it does not cause any human disease. It consists of a protein capsid (CAPside) and a 4.7 KB length single stranded DNA genome. The protein capsid consists of three subunits, VP1, VP2, and VP3.

Clinical and Translational Updates

Discounts, Gifts,
and more!

